



The Power of 3 =

*Faster speed + less fatigue*

Combining eye control, touch control, and EMG/Spatial control, the NeuroNode Trilogy is the latest system from Control Bionics providing a lifetime solution for those with paralysis and loss of speech. Faster communication, complete control powered by Grid 3, and a onetime funding application – what more could you ask for?

Say hello at [513-453-4848](tel:513-453-4848) and visit us at [controlbionics.com](http://controlbionics.com)

CONTROL BIONICS | LIFE CHANGING CONNECTIONS



#### WHAT IS THE NEURONODE TRILOGY?

The NeuroNode Trilogy is a 3-in-1 solution for those living with paralysis and loss of speech.

Combining touch control, eye control, and mind control, the NeuroNode Trilogy is finally the solution that fits to the users needs as they change. Combining control methods allows you to significantly increase your speed and communication output. In initial internal testing, users were able to communicate 133% faster with the NeuroNode Trilogy.

#### WHAT IS THE NEURONODE<sup>3</sup>?

At the heart of the NeuroNode Trilogy is the NeuroNode<sup>3</sup>, the latest generation of the world's leading wearable EMG control device. The NeuroNode<sup>3</sup> is packed with new added benefits – making it more intuitive and flexible than ever.

- Smaller, lighter, and even more durable
- Improved battery life for over 24 hours of use with hassle-free charging
- Adhesive-free, wearable solution with the new NeuroBand

The NeuroNode<sup>3</sup> now offers the combined power of its EMG + NEW 3D Spatial Control. This new control set allows users to choose EMG or Spatial control to access the NeuroNode Trilogy, all in one small, wearable device.

#### WHAT ARE MY FUNDING OPTIONS?

Our funding options are continuously expanding. Currently, we are proud to offer funding through Medicare, Medicaid, the Veterans Health Administration, as well as various Education Systems.

#### WHAT IS EMG?

EMG, or “electromyography,” is the measurement of electrical activity associated with the activation of a muscle group. While EMG technology has been used in clinical and research settings for decades, we have spent 10 years advancing the EMG technology inside the NeuroNode<sup>3</sup> specifically for people with severe motor impairments.

EMG control is a powerful access method for a wide range of conditions with many benefits:

- Entirely independent of body position and light condition
- Less tiring and more consistent than traditional access methods
- Completely wireless, making it easier for caregivers to work around

#### WHY CHOOSE THE NEURONODE TRILOGY?

While there are many great AAC solutions on the market, the NeuroNode Trilogy is designed to make your life easier. Combining three methods of control with Grid 3 software means one AAC software solution, one SLP evaluation, and one funding application. For super-users this means the fastest AAC solution around, and for those with progressive illnesses it means a solution that changes alongside their needs.

The Control Bionics Team supports the NeuroNode Trilogy around the clock. Give us a call any time of day and talk to a real person; most issues can be resolved within 3 minutes!

CONTROL BIONICS | LIFE CHANGING CONNECTIONS